

## SRRM1

---

**Reactivity:**Human Mouse Hamster

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:1000

**Calculated MW:**102kDa

**Observed MW:**Refer to Figures

**Immunogen:**

A synthetic peptide of human SRRM1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

160-KD; POP101; SRM160;

*To place an order, please [Click HERE](#).*

**Catalog #:**A6066

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**10250

**Isotype:**IgG

**Swiss Prot:**Q81YB3

**Purity:**Affinity purification

For research use only.