

## TIAL1

---

**Reactivity:** Human

**Tested applications:** WB

**Recommended Dilution:** WB 1:500 - 1:1000

**Calculated MW:** 41kDa

**Observed MW:** Refer to Figures

**Immunogen:**

Recombinant protein of human TIAL1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

j

**Synonym:**

TCBP; TIAR;

**Catalog #:** A6075

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 7073

**Isotype:** IgG

**Swiss Prot:** Q01085

**Purity:** Affinity purification

For research use only.

**Background:**

The protein encoded by this gene is a member of a family of RNA-binding proteins, has three RNA recognition motifs (RRMs), and binds adenine and uridine-rich elements in mRNA and pre-mRNAs of a wide range of genes. It regulates various activities including translational control, splicing and apoptosis. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. The different isoforms have been shown to function differently with respect to post-transcriptional silencing.

**To place an order, please [Click HERE](#).**