

ELAVL2

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:1000

Calculated MW:39kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of human ELAVL2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

HUB; HELN1; HEL-N1;

Catalog #:A6090

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:1993

Isotype:IgG

Swiss Prot:Q12926

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is a neural-specific RNA-binding protein that is known to bind to several 3' UTRs, including its own and also that of FOS and ID. The encoded protein may recognize a GAAA motif in the RNA. Three transcript variants encoding two different isoforms have been found for this gene.

To place an order, please [Click HERE](#).