

## KHDRBS2

---

**Reactivity:**Human

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:1000 IHC 1:100 - 1:200

**Calculated MW:**38kDa

**Observed MW:**Refer to Figures

**Immunogen:**

A synthetic peptide of human KHDRBS2

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

SLM1; SLM-1; bA535F17.1;

*To place an order, please [Click HERE](#).*

**Catalog #:**A6102

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**202559

**Isotype:**IgG

**Swiss Prot:**Q5VWX1

**Purity:**Affinity purification

For research use only.