

## HDAC11

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:**39kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human HDAC11

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

HD11;

**Catalog #:**A6140

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**79885

**Isotype:**IgG

**Swiss Prot:**Q96DB2

**Purity:**Affinity purification

For research use only.

**Background:**

This gene encodes a class IV histone deacetylase. The encoded protein is localized to the nucleus and may be involved in regulating the expression of interleukin 10. Alternative splicing results in multiple transcript variants.

*To place an order, please [Click HERE](#).*