

RAD54L2

Reactivity:Human

Tested applications:WB IF

Recommended Dilution:WB 1:500 - 1:2000 IF 1:20 - 1:100

Calculated MW:162kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human RAD54L2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ARIP4; HSPC325; SRISNF2L;

To place an order, please [Click HERE](#).

Catalog #:A6144

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:23132

Isotype:IgG

Swiss Prot:Q9Y4B4

Purity:Affinity purification

For research use only.