

VASH1

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 37kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human VASH1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

b

Synonym:

KIAA1036

Catalog #: A6148

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 22846

Isotype: IgG

Swiss Prot: Q7L8A9

Purity: Affinity purification

For research use only.

To place an order, please [Click HERE](#).