

DAP

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:20 - 1:100

Calculated MW:11kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human DAP

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #:A6150

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:1611

Isotype:IgG

Swiss Prot:P51397

Purity:Affinity purification

For research use only.

Background:

This gene encodes a basic, proline-rich, 15-kD protein. The protein acts as a positive mediator of programmed cell death that is induced by interferon-gamma. Alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.

To place an order, please [Click HERE](#).