

## HMGN2

**Reactivity:** Human Rat

**Tested applications:** WB IHC IF

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

**Calculated MW:** 9kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human HMGN2

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

f

**Synonym:**

HMG17;

**Catalog #:** A6156

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 3151

**Isotype:** IgG

**Swiss Prot:** P05204

**Purity:** Affinity purification

For research use only.

**Background:**

The protein encoded by this gene binds nucleosomal DNA and is associated with transcriptionally active chromatin. Along with a similar protein, HMGN1, the encoded protein may help maintain an open chromatin configuration around transcribable genes. The protein has also been found to have antimicrobial activity against bacteria, viruses and fungi.

*To place an order, please [Click HERE](#).*