

OSM

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 28kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human OSM

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A6163

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5008

Isotype: IgG

Swiss Prot: P13725

Purity: Affinity purification

For research use only.

Background:

Oncostatin M is a member of a cytokine family that includes leukemia-inhibitory factor, granulocyte colony-stimulating factor, and interleukin 6. This gene encodes a growth regulator which inhibits the proliferation of a number of tumor cell lines. It regulates cytokine production, including IL-6, G-CSF and GM-CSF from endothelial cells.

To place an order, please [Click HERE](#).