

## KLHL25

---

**Reactivity:**Human

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:2000

**Calculated MW:**66kDa

**Observed MW:**Refer to Figures

**Immunogen:**

Recombinant protein of human KLHL25

**Storage Buffer:**

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

ENC2; ENC-2;

*To place an order, please [Click HERE](#).*

**Catalog #:**A6167

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**64410

**Isotype:**IgG

**Swiss Prot:**Q9H0H3

**Purity:**Affinity purification

For research use only.