

SMYD5

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF IP

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:100 IF 1:50 - 1:100 IP 1:50 - 1:100

Calculated MW:47kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human SMYD5

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

f

Synonym:

RRG1; RAI15; NN8-4AG; ZMYND23

To place an order, please [Click HERE](#).

Catalog #:A6191

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:10322

Isotype:IgG

Swiss Prot:Q6GMV2

Purity:Affinity purification

For research use only.