

## TFPI2

**Reactivity:** Human

**Tested applications:** WB IHC

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:** 27kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human TFPI2

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

PP5; REF1; TFPI-2;

**Catalog #:** A6214

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 7980

**Isotype:** IgG

**Swiss Prot:** P48307

**Purity:** Affinity purification

For research use only.

**Background:**

This gene encodes a member of the Kunitz-type serine proteinase inhibitor family. The protein can inhibit a variety of serine proteases including factor VIIa/tissue factor, factor Xa, plasmin, trypsin, chymotrypsin and plasma kallikrein. This gene has been identified as a tumor suppressor gene in several types of cancer. Alternative splicing results in multiple transcript variants.

*To place an order, please [Click HERE](#).*