

ALPI

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 57kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ALPI

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

IAP;

Catalog #: A6226

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 248

Isotype: IgG

Swiss Prot: P09923

Purity: Affinity purification

For research use only.

Background:

There are at least four distinct but related alkaline phosphatases: intestinal, placental, placental-like, and liver/bone/kidney (tissue non-specific). The intestinal alkaline phosphatase gene encodes a digestive brush-border enzyme. This enzyme is a component of the gut mucosal defense system and is thought to function in the detoxification of lipopolysaccharide, and in the prevention of bacterial translocation in the gut.

To place an order, please [Click HERE](#).