

AKT1S1

Reactivity:Human Mouse

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:27kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human AKT1S1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

Lobe; PRAS40;

Background:

AKT1S1 is a proline-rich substrate of AKT (MIM 164730) that binds 14-3-3 protein (see YWHAH, MIM 113508) when phosphorylated (Kovacina et al., 2003)

To place an order, please [Click HERE](#).

Catalog #:A6238

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:84335

Isotype:IgG

Swiss Prot:Q96B36

Purity:Affinity purification

For research use only.