

## PNPLA2

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:**55kDa

**Observed MW:**Refer to Figures

**Immunogen:**

Recombinant protein of human PNPLA2

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

ATGL; TTS2; PEDF-R; FP17548; TTS-2.2; iPLA2zeta; 1110001C14Rik;

**Background:**

This gene encodes an enzyme which catalyzes the first step in the hydrolysis of triglycerides in adipose tissue. Mutations in this gene are associated with neutral lipid storage disease with myopathy.

*To place an order, please [Click HERE](#).*

**Catalog #:**A6245

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**57104

**Isotype:**IgG

**Swiss Prot:**Q96AD5

**Purity:**Affinity purification

For research use only.