

GAB1

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:20 - 1:100

Calculated MW:76kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human GAB1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #:A6248

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:2549

Isotype:IgG

Swiss Prot:Q13480

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the IRS1-like multisubstrate docking protein family. It is an important mediator of branching tubulogenesis and plays a central role in cellular growth response, transformation and apoptosis. Two transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).