

CDYL

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 66kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human CDYL

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CDYL1;

Catalog #: A6258

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9425

Isotype: IgG

Swiss Prot: Q9Y232

Purity: Affinity purification

For research use only.

Background:

Chromodomain Y is a primate-specific Y-chromosomal gene family expressed exclusively in the testis and implicated in infertility. Although the Y-linked genes are testis-specific, this autosomal gene is ubiquitously expressed. The Y-linked genes arose by retrotransposition of an mRNA from this gene, followed by amplification of the retroposed gene. Proteins encoded by this gene superfamily possess a chromodomain, a motif implicated in chromatin binding and gene suppression, and a catalytic domain believed to be involved in histone acetylation. Multiple proteins are encoded by transcript variants of this gene.

To place an order, please [Click HERE](#).