

TSPY3

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:1000 - 1:2000

Calculated MW:35kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of human TSPY3

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

f

Synonym:

CT78; TSPY1;TSPY;

To place an order, please [Click HERE](#).

Catalog #:A6270

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:728137

Isotype:IgG

Swiss Prot:P0CV98

Purity:Affinity purification

For research use only.