

## AGT

**Reactivity:** Human Mouse Rat

**Tested applications:** WB IHC

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:** 53kDa

**Observed MW:** Refer to Figures

**Immunogen:**

Recombinant protein of human AGT

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

ANHU; SERPINA8;

**Catalog #:** A6279

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 183

**Isotype:** IgG

**Swiss Prot:** P01019

**Purity:** Affinity purification

For research use only.

**Background:**

The protein encoded by this gene, pre-angiotensinogen or angiotensinogen precursor, is expressed in the liver and is cleaved by the enzyme renin in response to lowered blood pressure. The resulting product, angiotensin I, is then cleaved by angiotensin converting enzyme (ACE) to generate the physiologically active enzyme angiotensin II. The protein is involved in maintaining blood pressure and in the pathogenesis of essential hypertension and preeclampsia. Mutations in this gene are associated with susceptibility to essential hypertension, and can cause renal tubular dysgenesis, a severe disorder of renal tubular development. Defects in this gene have also been associated with non-familial structural atrial fibrillation, and inflammatory bowel disease.

**To place an order, please [Click HERE](#).**