

## BNIP2

---

**Reactivity:**Human Mouse

**Tested applications:**WB IHC IF

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

**Calculated MW:**36kDa

**Observed MW:**Refer to Figures

**Immunogen:**

Recombinant protein of human BNIP2

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

NIP2; BNIP-2;

**Background:**

This gene is a member of the BCL2/adenovirus E1B 19 kd-interacting protein (BNIP) family. It interacts with the E1B 19 kDa protein, which protects cells from virally-induced cell death. The encoded protein also interacts with E1B 19 kDa-like sequences of BCL2, another apoptotic protector.

*To place an order, please [Click HERE](#).*

**Catalog #:**A6282

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**663

**Isotype:**IgG

**Swiss Prot:**Q12982

**Purity:**Affinity purification

For research use only.