

PENK

Reactivity:Human Mouse

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW:30kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human PENK

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #:A6302

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:5179

Isotype:IgG

Swiss Prot:P01210

Purity:Affinity purification

For research use only.