

NOL3

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

Calculated MW:24kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human NOL3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ARC; FCM; MYP; NOP; NOP30;

Catalog #:A6319

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:8996

Isotype:IgG

Swiss Prot:O60936

Purity:Affinity purification

For research use only.

Background:

This gene encodes an anti-apoptotic protein that has been shown to down-regulate the enzyme activities of caspase 2, caspase 8 and tumor protein p53. Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).