

FAIM3

Reactivity:Human

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:43kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human FAIM3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

FCMR; TOSO;

Background:

Fc receptors specifically bind to the Fc region of immunoglobulins (Igs) to mediate the unique functions of each Ig class. FAIM3 encodes an Fc receptor for IgM (see MIM 147020) (Kubagawa et al., 2009

To place an order, please [Click HERE](#).

Catalog #:A6320

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:9214

Isotype:IgG

Swiss Prot:O60667

Purity:Affinity purification

For research use only.