

ZMYND11

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:70kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human ZMYND11

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

BS69; BRAM1;

Catalog #:A6327

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:10771

Isotype:IgG

Swiss Prot:Q15326

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene was first identified by its ability to bind the adenovirus E1A protein. The protein localizes to the nucleus. It functions as a transcriptional repressor, and expression of E1A inhibits this repression. Alternatively spliced transcript variants encoding different isoforms have been identified.

To place an order, please [Click HERE](#).