

CIAPIN1

Reactivity: Human Mouse

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 33kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human CIAPIN1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DRE2; PRO0915; Anamorsin; 2810413N20Rik;

Catalog #: A6336

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 57019

Isotype: IgG

Swiss Prot: Q6FI81

Purity: Affinity purification

For research use only.

Background:

CIAPIN1 is a cytokine-induced inhibitor of apoptosis with no relation to apoptosis regulatory molecules of the BCL2 (MIM 151430) or CASP (see MIM 147678) families. Expression of CIAPIN1 is dependent on growth factor stimulation (Shibayama et al., 2004)

To place an order, please [Click HERE](#).