

CANT1

Reactivity: Human

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100 IF 1:50 - 1:100

Calculated MW: 45kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CANT1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DBQD; SCAN1; SHAPY; SCAN-1;

Catalog #: A6341

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 124583

Isotype: IgG

Swiss Prot: Q8WVQ1

Purity: Affinity purification

For research use only.

Background:

This protein encoded by this gene belongs to the apyrase family. It functions as a calcium-dependent nucleotidase with a preference for UDP. Mutations in this gene are associated with Desbuquois dysplasia with hand anomalies. Alternatively spliced transcript variants have been noted for this gene.

To place an order, please [Click HERE](#).