

## STX6

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC IF IP

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:20 - 1:200 IF 1:20 - 1:100 IP 1:20 - 1:50

**Calculated MW:**29kDa

**Observed MW:**Refer to figures

**Immunogen:**

A synthetic Peptide of human STX6

**Storage Buffer:**

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

*To place an order, please [Click HERE](#).*

**Catalog #:**A6343

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**10228

**Isotype:**IgG

**Swiss Prot:**O43752

**Purity:**Affinity purification

For research use only.