

PDLIM3

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:39kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human PDLIM3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ALP;

Catalog #:A6346

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:27295

Isotype:IgG

Swiss Prot:Q53GG5

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene contains a PDZ domain and a LIM domain, indicating that it may be involved in cytoskeletal assembly. In support of this, the encoded protein has been shown to bind the spectrin-like repeats of alpha-actinin-2 and to colocalize with alpha-actinin-2 at the Z lines of skeletal muscle. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. Aberrant alternative splicing of this gene may play a role in myotonic dystrophy.

To place an order, please [Click HERE](#).