

## BTG1

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC IF

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

**Calculated MW:**19kDa

**Observed MW:**Refer to Figures

**Immunogen:**

Recombinant protein of human BTG1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Catalog #:**A6359

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**694

**Isotype:**IgG

**Swiss Prot:**P62324

**Purity:**Affinity purification

For research use only.

**Background:**

This gene is a member of an anti-proliferative gene family that regulates cell growth and differentiation. Expression of this gene is highest in the G0/G1 phases of the cell cycle and downregulated when cells progressed through G1. The encoded protein interacts with several nuclear receptors, and functions as a coactivator of cell differentiation. This locus has been shown to be involved in a t(8;12)(q24;q22) chromosomal translocation in a case of B-cell chronic lymphocytic leukemia.

*To place an order, please [Click HERE](#).*