

GP2

Reactivity:Human Mouse

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW:59kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human GP2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ZAP75;

Catalog #:A6378

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:2813

Isotype:IgG

Swiss Prot:P55259

Purity:Affinity purification

For research use only.

Background:

This gene encodes an integral membrane protein that is secreted from intracellular zymogen granules and associates with the plasma membrane via glycosylphosphatidylinositol (GPI) linkage. The encoded protein binds pathogens such as enterobacteria, thereby playing an important role in the innate immune response. The C-terminus of this protein is related to the C-terminus of the protein encoded by the neighboring gene, uromodulin (UMOD). Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).