

## PNLIP

**Reactivity:** Human Mouse Rat

**Tested applications:** WB IHC IF

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

**Calculated MW:** 51kDa

**Observed MW:** Refer to Figures

**Immunogen:**

Recombinant protein of human PNLIP

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

PL; PTL; PNLIPD;

**Catalog #:** A6396

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 5406

**Isotype:** IgG

**Swiss Prot:** P16233

**Purity:** Affinity purification

For research use only.

**Background:**

This gene is a member of the lipase gene family. It encodes a carboxyl esterase that hydrolyzes insoluble, emulsified triglycerides, and is essential for the efficient digestion of dietary fats. This gene is expressed specifically in the pancreas.

*To place an order, please [Click HERE](#).*