

PPOX

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW:50kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human PPOX

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

VP; PPO; V290M;

Background:

This gene encodes the penultimate enzyme of heme biosynthesis, which catalyzes the 6-electron oxidation of protoporphyrinogen IX to form protoporphyrin IX. Mutations in this gene cause variegate porphyria, an autosomal dominant disorder of heme metabolism resulting from a deficiency in protoporphyrinogen oxidase, an enzyme located on the inner mitochondrial membrane. Alternatively spliced transcript variants encoding the same protein have been identified.

To place an order, please [Click HERE](#).

Catalog #:A6397

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:5498

Isotype:IgG

Swiss Prot:P50336

Purity:Affinity purification

For research use only.