

RCVRN

Reactivity:Human

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

Calculated MW:23kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human RCVRN

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RCV1;

Catalog #:A6404

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:5957

Isotype:IgG

Swiss Prot:P35243

Purity:Affinity purification

For research use only.

Background:

This gene encodes a member of the recoverin family of neuronal calcium sensors. The encoded protein contains three calcium-binding EF-hand domains and may prolong the termination of the phototransduction cascade in the retina by blocking the phosphorylation of photo-activated rhodopsin. Recoverin may be the antigen responsible for cancer-associated retinopathy.

To place an order, please [Click HERE](#).