

## SFRP4

**Reactivity:** Human Mouse Rat

**Tested applications:** WB IHC IF

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:100 IF 1:50 - 1:100

**Calculated MW:** 39kDa

**Observed MW:** Refer to Figures

**Immunogen:**

Recombinant protein of human SFRP4

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

FRP-4; FRPHE; sFRP-4;

**Catalog #:** A6409

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 6424

**Isotype:** IgG

**Swiss Prot:** Q6FHJ7

**Purity:** Affinity purification

For research use only.

**Background:**

Secreted frizzled-related protein 4 (SFRP4) is a member of the SFRP family that contains a cysteine-rich domain homologous to the putative Wnt-binding site of Frizzled proteins. SFRPs act as soluble modulators of Wnt signaling. The expression of SFRP4 in ventricular myocardium correlates with apoptosis related gene expression.

*To place an order, please [Click HERE](#).*