

SFRP4

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:100 IF 1:50 - 1:100

Calculated MW:39kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human SFRP4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

FRP-4; FRPHE; sFRP-4;

Catalog #:A6409

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:6424

Isotype:IgG

Swiss Prot:Q6FHJ7

Purity:Affinity purification

For research use only.

Background:

Secreted frizzled-related protein 4 (SFRP4) is a member of the SFRP family that contains a cysteine-rich domain homologous to the putative Wnt-binding site of Frizzled proteins. SFRPs act as soluble modulators of Wnt signaling. The expression of SFRP4 in ventricular myocardium correlates with apoptosis related gene expression.

To place an order, please [Click HERE](#).