

## TCN1

---

**Reactivity:**Human

**Tested applications:**WB IHC IF

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

**Calculated MW:**48kDa

**Observed MW:**Refer to Figures

**Immunogen:**

Recombinant protein of human TCN1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

HC; TC1; TCI; TC-1;

**Background:**

This gene encodes a member of the vitamin B12-binding protein family. This family of proteins, alternatively referred to as R binders, is expressed in various tissues and secretions. This protein is a major constituent of secondary granules in neutrophils and facilitates the transport of cobalamin into cells.

*To place an order, please [Click HERE](#).*

**Catalog #:**A6414

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**6947

**Isotype:**IgG

**Swiss Prot:**P20061

**Purity:**Affinity purification

For research use only.