

## SF1

**Reactivity:** Human Mouse Rat

**Tested applications:** WB IHC

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:** 68kDa

**Observed MW:** Refer to Figures

**Immunogen:**

Recombinant protein of human SF1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

BBP; MBBP; ZFM1; ZNF162; D11S636; ZCCHC25;

**Catalog #:** A6424

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 7536

**Isotype:** IgG

**Swiss Prot:** Q15637

**Purity:** Affinity purification

For research use only.

**Background:**

This gene encodes a nuclear pre-mRNA splicing factor. The encoded protein specifically recognizes the intron branch point sequence and is required for the early stages of spliceosome assembly. Alternate splicing results in multiple transcript variants.

*To place an order, please [Click HERE](#).*