

FIBP

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 41kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human FIBP

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

FGFIBP; FIBP-1;

Catalog #: A6436

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9158

Isotype: IgG

Swiss Prot: O43427

Purity: Affinity purification

For research use only.

Background:

Acidic fibroblast growth factor is mitogenic for a variety of different cell types and acts by stimulating mitogenesis or inducing morphological changes and differentiation. The FIBP protein is an intracellular protein that binds selectively to acidic fibroblast growth factor (aFGF). It is postulated that FIBP may be involved in the mitogenic action of aFGF. Two transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).