

C1D

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW:16kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human C1D

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

LRP1; hC1D; SUNCOR; SUN-CoR;

Catalog #:A6448

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:10438

Isotype:IgG

Swiss Prot:Q13901

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is a DNA binding and apoptosis-inducing protein and is localized in the nucleus. It is also a Rac3-interacting protein which acts as a corepressor for the thyroid hormone receptor. This protein is thought to regulate TRAX/Translin complex formation. Alternate splicing results in multiple transcript variants that encode the same protein. Multiple pseudogenes of this gene are found on chromosome 10.

To place an order, please [Click HERE](#).