

LYVE1

Reactivity:Human Mouse

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:35kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human LYVE1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HAR; XLKD1; LYVE-1; CRSBP-1;

Background:

This gene encodes a type I integral membrane glycoprotein. The encoded protein acts as a receptor and binds to both soluble and immobilized hyaluronan. This protein may function in lymphatic hyaluronan transport and have a role in tumor metastasis.

To place an order, please [Click HERE](#).

Catalog #:A6452

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:10894

Isotype:IgG

Swiss Prot:Q9Y5Y7

Purity:Affinity purification

For research use only.