

SAP30BP

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 33kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human SAP30BP

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HTRG; HTRP; HCNGP;

Catalog #: A6460

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 29115

Isotype: IgG

Swiss Prot: Q9UHR5

Purity: Affinity purification

For research use only.

Background:

SAP30BP (SAP30 binding protein) is a protein-coding gene. Diseases associated with SAP30BP include obesity, and schizophrenia.

To place an order, please [Click HERE](#).