

CPA6

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:100 - 1:200

Calculated MW: 51kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CPA6

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CPAH; ETL5; FEB11;

Catalog #: A6475

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 57094

Isotype: IgG

Swiss Prot: Q8N4T0

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene belongs to the family of carboxypeptidases, which catalyze the release of C-terminal amino acid, and have functions ranging from digestion of food to selective biosynthesis of neuroendocrine peptides. Polymorphic variants and a reciprocal translocation t(6;8)(q26;q13) involving this gene, have been associated with Duane retraction syndrome.

To place an order, please [Click HERE](#).