

## UCN2

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC IF

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:100 IF 1:50 - 1:100

**Calculated MW:**12kDa

**Observed MW:**Refer to Figures

**Immunogen:**

Recombinant protein of human UCN2

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

UR; SRP; URP; UCNI; UCN-II;

**Catalog #:**A6485

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**90226

**Isotype:**IgG

**Swiss Prot:**Q96RP3

**Purity:**Affinity purification

For research use only.

**Background:**

This gene is a member of the sauvagine/corticotropin-releasing factor/urotensin I family. It is structurally related to the corticotropin-releasing factor (CRF) gene and the encoded product is an endogenous ligand for CRF type 2 receptors. In the brain it may be responsible for the effects of stress on appetite. In spite of the gene family name similarity, the product of this gene has no sequence similarity to urotensin II.

*To place an order, please [Click HERE](#).*