

GRK1

Reactivity:Human Mouse

Tested applications:WB IHC

Recommended Dilution:WB 1:200 - 1:2000 IHC 1:50 - 1:200

Calculated MW:64kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of human GRK1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RK; RHOK; GPRK1;

Catalog #:A6497

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:6011

Isotype:IgG

Swiss Prot:Q15835

Purity:Affinity purification

For research use only.

Background:

This gene encodes a member of the guanine nucleotide-binding protein (G protein)-coupled receptor kinase subfamily of the Ser/Thr protein kinase family. The protein phosphorylates rhodopsin and initiates its deactivation. Defects in GRK1 are known to cause Oguchi disease 2 (also known as stationary night blindness Oguchi type-2).

To place an order, please [Click HERE](#).