

ARL3

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:20kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic Peptide of human ARL3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

1 mg/ml

Synonym:

ARFL3;

Catalog #:A6499

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:403

Isotype:IgG

Swiss Prot:P36405

Purity:Affinity purification

For research use only.

Background:

ADP-ribosylation factor-like 3 is a member of the ADP-ribosylation factor family of GTP-binding proteins. ARL3 binds guanine nucleotides but lacks ADP-ribosylation factor activity.

To place an order, please [Click HERE](#).