

PRDM14

Reactivity:Mouse

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:64kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human PRDM14

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

To place an order, please [Click HERE](#).

Catalog #:A6502

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:63978

Isotype:IgG

Swiss Prot:Q9GZV8

Purity:Affinity purification

For research use only.