

ASF1B

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:200 - 1:1000

Calculated MW:22kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human ASF1B

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CIA-II;

Catalog #:A6511

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:55723

Isotype:IgG

Swiss Prot:Q9NVP2

Purity:Affinity purification

For research use only.

Background:

This gene encodes a member of the H3/H4 family of histone chaperone proteins and is similar to the anti-silencing function-1 gene in yeast. The encoded protein is the substrate of the tousel-like kinase family of cell cycle-regulated kinases, and may play a key role in modulating the nucleosome structure of chromatin by ensuring a constant supply of histones at sites of nucleosome assembly.

To place an order, please [Click HERE](#).