

ADSS

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

Calculated MW:50kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human ADSS

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ADEH; ADSS 2;

Catalog #:A6516

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:159

Isotype:IgG

Swiss Prot:P30520

Purity:Affinity purification

For research use only.

Background:

This gene encodes the enzyme adenylosuccinate synthetase which catalyzes the first committed step in the conversion of inosine monophosphate to adenosine monophosphate. A pseudogene of this gene is found on chromosome 17.

To place an order, please [Click HERE](#).