

## ATOH1

---

**Reactivity:**Human

**Tested applications:**WB IHC IF

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

**Calculated MW:**38kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human ATOH1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

ATH1; HATH1; MATH-1; bHLHa14;

**Background:**

This protein belongs to the basic helix-loop-helix (BHLH) family of transcription factors. It activates E-box dependent transcription along with E47.

*To place an order, please [Click HERE](#).*

**Catalog #:**A6530

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**474

**Isotype:**IgG

**Swiss Prot:**Q92858

**Purity:**Affinity purification

For research use only.